

The **CS202I-DMX-4SS** provides the same high performance as our CS202I-DMX-1SS with its all welded stainless steel construction and welded stainless steel instrumentation skirt but provides a larger than standard sample space. The vacuum shroud comes standard with 4 window ports, however a 5th port can be added on the end. The system is capable of vacuum levels of 10<sup>-7</sup> Torr with an appropriate vacuum pump.

# **Applications**

- Large Samples
- Optical
- Raman
- UV, VIS, IR
- FTIR
- Electro & Photoluminescence
- Resistivity/Hall Probe Experiments
- Diamond Anvil Cell
- Magneto-Optical
- PITS / DLTS
- Thermal, Electrical and Magnetic Susceptibility
- Magneto Optical Kerr Effect (MOKE)

# Features

- Cryogen Free, Low Power
- High Performance Stainless Steel Construction
- Large clear view optical windows (1.5 in)
- Large sample viewing angle for optical collection (F/1.1)
- Can operate in any orientation
- Fully customizable

# **Typical Configuration**

- Cold head (DE-202AI)
- Compressor (ARS-4HW)
- 2 Helium Hoses
- Stainless Steel vacuum shroud with 4 window ports for optical and electrical measures (DMX-4SS)
- Nickel Plated OFHC radiation shield
- 2 High purity quartz windows
- Instrumentation for temperature measurement and control:
  - 10 pin hermetic feed through 36 ohm thermofoil heater Silicon diode sensor curve matched to (±0.5K) for control Calibrated silicon diode sensor (±12 mk) with 4 in. free length for accurate sample measurement.
- Wiring for electrical experiments: 10 pin hermetic feed through 4 copper wires
- Sample holder for optical and electrical experiments
- Temperature Controller

# **Options and Upgrades**

- 4K Coldhead (0.1W @ 4.2K)
- 5.5K Coldhead (1W @ 10K)
- 450K High Temperature Interface
- 800K High Temperature Interface
- Turbo upgrade for faster cooldown times
- Custom temperature sensor configuration (please contact our sales staff
- Custom wiring configurations (please contact our sales staff)
- Window material upgrades (custom materials available)
- Sample holder upgrades (custom sample holders available)



The above picture shows a cryocooler with a vacuum shroud, radiation shield, and sample holder installed.



The above picture shows a complete system (minus the vacuum pump and temperature controller).

# www.arscryo.com



# **Optical Cryostat - Large Sample Space**

# **Cooling Technology**

	DE-202	Closed Cycle Cryocooler		
	Refrigeration Type	Pneumatically Driven GM Cycle		
	Liquid Cryogen Usage	None, Cryogen Free		
Tem	iperature*			
	DE-202AI	< 10K - 350K		
	DE-202PI	< 5.5K - 350K		
	DE-202SI	< 4.2K - 350K		
	With 800K Interface	(Base Temp + 2K) - 700K		
	With 450K Interface	(Base Temp + 2K) - 450K		
	Stability	0.1K		

 $^{*}\mbox{Based}$  on bare cold head with a closed radiation shield, and no additional sources of experimental or parasitic heat load

#### Sample Space

Ор

	Diameter	47.5 mm (1.87in.)				
	Height	90 mm (3.54 in.)				
	Sample Holder Attachment	1/4 - 28 screw				
	Sample Holder	www.arscryo.com/Products/ SampleHolders.html				
oti	ical Access					
	Window Ports	4 - 90° Apart				
	Diameter	51 mm (2 in)				

	Diameter	51 mm (2 in)				
	Clear View	38 mm (1.5 in)				
#/F		1.1				
Window Material		www.arscryo.com/Products/ WindowMaterials.html				


### Temperature Instrumentation and Control (Standard)

Heater		36 ohm Thermofoil Heater anchored to the coldtip					
	Control Sensor	Curve Matched Silicon Diode installed on the coldtip					
	Sample Sensor	Calibrated Silicon Diode with free length wires					
	Contact ARS for other options						
Inst							
	Instrumentation Skirt	Welded Stainless Steel					
	Pump out Port	1 - NW 25					
	Instrumentation Ports	3					
	Instrumentation Wiring	Contact sales staff for options					
Vac	uum Shroud						
	Material	Stainless Steel					
	Length	378 mm (14.9 in)					
	Diameter	95mm (3.75 in) (at the sample space)					
	Width	95 mm (3.75 in) (at the sample space)					
Rad	iation Shield						
	Material	Nickel Plated OFHC Copper					
	Attachment	Flanged					
	Optical Access	0, 2, or 4 (customer specified)					
Cryostat Footprint							
	Overall Length	603 mm (23.72 in)					
	Motor Housing Diameter	114 mm (4.5 in)					
	Rotational Clearance	200 mm (8 in) with "G" Configuration					

Cryocooler Model		DE-2	DE-202AI		DE-202A(T)I		DE-202PI		DE-202SI	
	Frequency	60 Hz	50 Hz	60 Hz	50 Hz	60 Hz	50 Hz	60 Hz	50 Hz	
Base Temperature		<9K	<9K	<9K	<9K	<5.5K	<5.5K	<4.2K	<4.2K	
Cooling Capacity*	4.2K	-	-	-	-	-	-	0.1W	0.08W	
	10K	0.5W	0.4W	0.7W	0.56W	1W	0.8W	1.2W	1W	
	20К	2.5W	2W	3.7W	3W	3.5W	2.8W	4W	3.2W	
	77K	4W	3.2W	6W	4.8W	3.5W	2.8W	4W	3.2W	
Radiation Shield Co	ooling Capacity	10W	8W	15W	12W	10W	8W	10W	8W	
Cooldown Time	20K	50 min	60 min	35 min	42 min	60 min	72 min	60 min	72 min	
	Base Temperature	70 min	84 min	50 min	60 min	90 min	108 min	90 min	108 min	
Compressor Model		ARS-	4HW	ARS-	4HW	ARS-	4HW	ARS-	4HW	
Typical Maintenand	ce Cycle	12,000	) hours	12,000	hours	12,000	) hours	12,000	hours	



# CS202\*I-DMX-4SS Outline Drawing







# Optical Cryostat - Large Sample Space

# Sample Space

### **ARS-4HW Compressor**





Compresso	r Model	ARS-4HW			
	Frequency	60 Hz	50 Hz		
Standard Voltage	Min	208 V	190 V		
	Max	230 V	210 V		
Transformer Options	10%		220 V, 230 V		
	15%		240 V		
Power Usage	Single Phase	3.6 kW	3.0 kW		
Refrigerant Gas		99.999% Helium Gas, Pre-Charged			
Noise Level		60 dBA			
Ambient Temperature		12 - 40 C (54 - 104 F)			
Cooling Water Consumption		2.3 L / min (0.6 Gal. / min)			
Temperature		10 - 35 C (50–95 F)			
Connection		3/8 in. Swagelok Fitting			
Dimensions:	L	483 mm (19 in)			
W		434 mm (17.1 in)			
	Н	516 mm (20.3 in)			
Weight		72 kg (160 lbs)			
Typical Maintenance Cyc	le	12,000 hours			
Water Recirculation Opt	ion	CoolPac Compatible			